

SEATING



The Pedal Rocker is a superior mobility chair that provides a gentle rocking motion.

STANDARD FEATURES

- Broda Comfort Tension Seating™
- 16 gauge tubular steel powder coated frame
- Rocking function with caregiver lock
- Adjustable rocking tension to accommodate various residents
- Up to 18° of infinitely adjustable seat tilt (Deep tilt position of 33° can be achieved when tilt is combined with rocking function)
- · Adjustable arm height
- Removable padding including neckrest and lower side pads
- 18" or 20" standard seat width
- 14" to 19" adjustable seat height
- 17" seat depth
- 32" back length
- 5" total lock front casters, 5" fixed rear casters
- 250 lb. weight capacity

PEDAL ROCKER

WHERE COMFORT IS ESSENTIAL

- Residents can self propel and rock in the same chair! Excellent for those who are mobile that may be agitated and have a tendency to rock.
- Rocking motion helps to calm individuals and helps prevent wandering Excellent for Alzheimers. Rocking also helps absorb rigid movements to aid in resident safety and comfort.
- 18° of infinitely adjustable tilt in space cradles the body helping to prevent slumping, sliding and falls without the use of restraints. It also provides exceptional comfort while helping to reduce pressure.
- Front pivot seat tilt allows for a proper foot on floor position for foot propelling while tilted Encourages safe mobility and independence.
- Adjustable seat to floor height from 14" to 19" accommodates both short and tall residents.



Pedal Rocker with full padding package



Resident can use rocking function independently of tilt (Rocking function shown)

NEW! PEDAL ROCKER

EXCELLENT FOR FALL PREVENTION

OPTIONS

- Full padding package including neckrest, wings, shoulder bolsters, lower side pads, and seat and back with removable terry or vinyl covers. Full vinyl covers also available (right) as additional option.
- 20" or 22" mag wheels
- Solid or clear tray
- Padded tray cover
- IV Pole
- Molded laterals
- 15.5", 18.5", 20" seat depth
- Swing away and elevating swing away footrests
- Half tray
- 02 Holder
- Thigh belt
- 16" or 22" seat width
- 33.5", 35", 36", 38", 40" back length

The ergonomic design maintains a proper foot on floor position when using the seat tilt, helping to prevent falls.



18° of seat tilt



Full seat tilt combined with rocking 33°





SEATING



When not self propelling, residents can use the tilt function along with the rocking function to achieve a deep seat tilt position. Excellent for napping!

Broda's Comfort Tension Seating™

The Key to Comfort is in Broda's Strapping! Memory retention allows straps to return to original shape within seconds.

- Redistributes pressure helping to prevent skin breakdown
- Creates air flow for heat and moisture exchange
- Contours to the body for long term sitting comfort

