

*Seating is Believing*

*Providing superior quality and comfort for a variety of needs*



• **Pressure Redistribution**

• **Long Term Comfort**

• **Restraint Reduction**

• **Fall Prevention**

### **Broda's Comfort Tension Seating™ is the Key to Comfort!**

- Redistributes pressure to help prevent skin breakdown
- Creates air flow for heat and moisture exchange
- Contours to the body for long term sitting comfort
- Memory retention allows straps to return to their original shape within seconds for long term and multi resident use

**Broda Comfort Tension Seating™**



Since 1981, Broda Seating has been manufacturing quality repositioning chairs. Our chairs help improve quality of life as well as help ease caregiver burden. Known for our comfort, pressure redistribution and durability, Broda chairs are an integral part of any facility in that we:

- Help to prevent skin breakdown
- Aid in Safe Patient Handling and No-Fall programs, thus helping to reduce injuries
- Accommodate difficult to seat residents
- Help to reduce chair maintenance and replacement costs



## POSITIONING CHAIRS

- Excellent restraint-free, long term positioning
- Suit a variety of needs from general comfort, to difficult positioning issues.



30Vt with vinyl seat and back pad upgrade

### ▲ Basic Chairs Models 30VT & 80V

- Basic Chairs are ideal for hospice care and general resident seating.
- An excellent cost effective alternative to the Midline and Elite models
- Both 30VT and 80V offer seat tilt of 25°, Model 30VT offers additional positioning with 45° of back recline
- Standard Basic Padding Package: Neckrest, Shoulder Bolsters, Lower Side Pads, Sole Pad



785 Tilt Recliner

### ▶ Midline Full Recliner

- The Midline offers a full range of positioning on a budget, with 38° of seat tilt and 90° of full back recline.
- Adjustable seat height of 17"–22" allows for taller residents to use as a mobility chair with the legrest removed.
- Standard Full Padding Package: as shown



Midline Full Recliner

### ▼ Elite Chairs Models 785 and 85V

- Elite Chairs offer a full range of positioning and function, and are ideal for a wide range of needs including your most difficult to seat residents.
- Both 785 and 85V offer seat tilt of 38°. Model 785 offers additional positioning with 90° of full back recline
- Swing away removable arms and wings assist with transfers.
- 6" multi-directional rear casters offer easy one-step function
- Ideal for severe agitation and involuntary movement. Huntington's Padding (HSP Padding Package) offers additional safety and comfort.
- Standard Full Padding Package: as shown

## GLIDERS & MOBILITY CHAIRS

- Encourage mobility and independence
- Pedal models offer front pivot seat tilt for a proper "foot on floor" position for safe propelling



100-10AL Adult Small Glider with optional wooden arms

### ◀ Auto Locking Glider Models Tall Back 100-20AL, Low Back 100-15AL, Adult Small 100-10AL

- The Auto Locking Glider is excellent for resident rooms and common areas
- Ideal for Alzheimer's as it produces a calming effect
- Gentle gliding helps reduce the need to wander
- Chair begins to glide once user sits down and stops gliding once user begins to stand



*Pedal Chair with full padding upgrade*

## ◀ Pedal Chair *Model 48*

- The Pedal Chair is the ideal alternative to a wheelchair
- 18° of seat tilt helps prevent falls and allows for safe propelling
- Adjustable seat height of 12.5" -17.5" accommodates short or tall residents
- Standard Pedal Padding: Neckrest, Lower Side pads

## Pedal Rocker *Model 48R* ▶

- The Pedal Rocker offers both mobility and gentle rocking
- Rocking action is ideal for Alzheimer's
- Up to 33° seat tilt when tilt is combined with rocking function
- Adjustable seat height of 14"-19"
- Standard Pedal Padding: Neckrest, Lower Side Pads



*Pedal Rocker with full padding upgrade*

- See Bariatric Chairs below for our Bariatric Commode Shower Chair

## SHOWER COMMODE

## ◀ CS 385 Commode Shower

- The CS 385 offers superior function and comfort in a commode shower chair
- 12° of posterior seat tilt, and 10° of anterior tilt help prevent falls and ease caregiver burden
- 67° of back recline allows for easy cleaning and washing resident's hair
- Soft commode seat provides exceptional comfort



*CS 385 with optional headrest and laterals*

## BARIATRIC CHAIRS

- Superior durability, positioning and comfort for bariatric care
- Excellent for Safe Patient Handling

## ◀ Bari 385 Commode Shower *Model B385*

- The Bari 385 offers a full range of function for bariatric commode and shower needs
- 800 lb. weight capacity
- Adjustable seat width of 27" - 33"
- 14° of posterior tilt and 8° of anterior tilt
- 65° of back recline helps with washing hard to reach areas



*Bari 385 with optional headrest*

## Vanguard Elite Bariatric *Model 985* ▶

- The Vanguard offers outstanding comfort unlike any other bariatric chair
- 500 lb. weight capacity
- 30° of seat tilt and 10° anterior tilt operated by foot pump- for No Lift
- 45° of back recline assists with comfort and breathing
- 26", 28", 30" seat widths



*Vanguard Elite with IV pole and vinyl seat and back pad upgrade.*

**Broda chairs feature infinitely adjustable tilt-in-space to allow for an infinite range of positioning options for optimal comfort, safety and repositioning.**




# PRODUCT SPECIFICATIONS

|  | POSITIONING CHAIRS  | MOBILITY CHAIRS  | GLIDERS  | SHOWER COMMODE  | BARIATRIC CHAIRS  |
|--|---|--|--|---|---|
| <b>Tilt (Infinitely Adjustable)</b>    | <b>Elite Chairs 785 and 85V</b><br>Model 785 up to 38°<br>Model 85V up to 38°<br>Model 785 up to 90° (full recline)<br>Model 85V up to 45° (Semi recline 45° available)<br>Model 85V no recline | <b>48 Pedal Chair</b><br>Up to 18°<br>No Recline                                   | <b>Auto Locking Glider Models 100-20AL, 100-15AL, 100-10AL (Adult Small)</b><br>N/A<br>N/A | <b>CS 385</b><br>Up to 12° posterior<br>Up to 10° anterior<br>Up to 67° | <b>Bari 385</b><br>Up to 14° posterior<br>Up to 8° anterior<br>Up to 65°  |
| <b>Recline (Infinitely Adjustable)</b> | Model 30VT up to 25°<br>Model 80V up to 25°<br>Model 30VT up to 45°<br>Model 80V no recline   | <b>48R Pedal Rocker</b><br>Up to 18°<br>Up to 33° total when combined with rocking | N/A  | Up to 12° posterior<br>Up to 10° anterior                               | Vanguard Elite 985 Tilt Recliner<br>Up to 30° posterior<br>Up to 10° anterior/tilt                              |
| <b>Frame</b>                           | 16 gauge tubular steel  | 16 gauge tubular steel   | 16 gauge tubular steel   | 16 gauge stainless steel  | 16 gauge stainless steel  |
| <b>Standard Padding</b>                | Neckrest, Shoulder Bolsters, Upper Side Pads, Lower Side Pads, Sole Pad   | Neckrest, Lower Side Pads  | Neckrest, Lower Side Pads, Seat and Back Pads  | N/A   | Neckrest, Shoulder Bolsters, Upper Side Pads, Lower Side Pads, Legrest Pad, Seat and Back Pad with Terry Covers |
| <b>Standard Seat Width</b>             | 18" or 20"  | 18" or 20"   | 20" (100-20AL, 100-15AL) 18" (100-10AL)  | 22"   | 27" (33" with arm width adjustment)   |
| <b>Standard Seat Depth</b>             | 17"   | 17"  | 17"  | 18"   | Adjustable 21" or 23"<br>Concave back for gluteal shelf   |
| <b>Standard Seat to Floor</b>          | Adjustable 19.5", 20.5", 21.5"  | Adjustable from 12.5" to 17.5"   | 18" (100-20AL, 100-15AL) 16" (100-10AL)  | 24.5" (rolls over most standard toilets)                                | Adjustable 21" or 23"<br>Optional Concave back for gluteal shelf  |
| <b>Standard Back Length</b>            | 32"   | 32"  | 32" (100-20AL) 26" (100-15AL, 100-10AL)  | 26"   | Adjustable 20", 21", 22", 23", 24"<br>Adjustable 32", 34", 36", 38", 40"  |
| <b>Custom Seat Widths</b>              | 16", 22", 24", 26", 28"   | 16", 22", 24", 26", 28"  | N/A  | Optional extended arms achieve 25" seat width                           | N/A   |
| <b>Arms/Side Panels</b>                | Height adjustable with push pins, swing away & removable side panels  | Height adjustable with push pins   | Height Adjustable with push pins   | Swing away and removable  | Removable, adjustable fore and aft positions  |
| <b>Legrest</b>                         | Elevates with back recline (and extends with recline on 785)<br>Elevates independently Length adjustable  | N/A  | N/A  | Elevates independently of back recline Length adjustable                | Elevates independently of back recline, and independently Length Adjustable                                     |
| <b>Footrest</b>                        | Split or Single Flip Down   | Upgrade only: Swing Away or Elevating Swing Away                                   | N/A  | Single Flip Down  | Single Flip Down  |
| <b>Casters</b>                         | 2 x 5" Front Total Lock<br>2 x 6" Rear Multi-Directional  | 2 x 5" Front Total Lock<br>2 x 5" Rear Fixed                                       | 2 rear transport wheels  | 3 x 5" Total Lock Stainless<br>1 x 5" Directional Lock Stainless        | 2 x 6" Front Total Lock<br>2 x 6" Rear Multi-Directional  |
| <b>Weight Capacity</b>                 | 350 lbs.  | 350 lbs.   | 250 lbs.   | 350 lbs.  | 800 lbs.<br>500 lbs.  |
| <b>Options</b>                         | Yes   | No   | No   | No  | No  |
| <b>HSP (Huntington's) Padding*</b>     | Yes   | No   | No   | No  | No  |
| <b>Mag Wheels</b>                      | Yes- rear mounted only  | Yes- front or rear mounted   | No   | No  | No  |
| <b>Contracture Sling</b>               | Model 785 - Yes- semi recline only<br>Model 85V - Yes   | Yes  | No   | No  | No  |

## See Broda Accessories List for additional options

HSP \*Short Wing Covers, Extended Shoulder Bolsters, Thigh Belt, Full Extended Footpad, Fixed Footrest (factory installed)



Shown with full padding upgrade

Shown with full padding upgrade

Shown with optional headrest

1-800-668-0637